

SERUM STICK

Transformation begins within,
Tap tap - reveal your radiance, anytime, anywhere

Revitalize Your Skin with Intensive,
High-Concentration Serums.



SERUM STICK

Easy to carry, fits right in your pocket!

Choose the serum you want and apply it instantly to your skin!
The perfect daily care tailored for you!

SERUM STICK

Delivering serum deep into the skin with **EP technology**
— simple, powerful care in a **convenient stick!**

The Serum Stick offers a convenient way to apply highly concentrated serum in a stick form. With EP technology, it maximizes skin absorption, providing deep and powerful care.

Its specialized serums improve skin tone and elasticity, delivering moisture and nutrients effectively. Suitable for all skin types, it's hygienic, easy to use, and delivers professional-quality skincare at home.



Patents Pending





Skincare Device Equipped with a Medication Container (Application No. 10-2025-0031654)
Electroporation-Based Drug Delivery Skincare Device (Application No. 10-2025-0031932)

EP (Electroporation) Electroporation technology

An innovative skincare technology that uses microcurrents to temporarily open the skin's natural barrier, allowing highly concentrated serums and active ingredients to be absorbed deep within the skin.

Precise 3-level adjustment

Customizable based on individual skin condition and sensitivity




Just one Serum Stick!

ALL AT ONCE!!


Up to **500x**  Greater Absorption Power




Absorption Boost with EP



VS
<



VS
≪



Compared to hands or galvanic devices, the Serum Stick shows significantly higher absorption.

* According to an SCI-level study published in the Journal of Investigative Dermatology (1998), EP (Electroporation) enhances absorption efficiency up to 500 times compared to manual application.


Tested using a Raman Spectrophotometer

Skin Absorption Test (absorption amount, depth, and speed)

Before application

0 μ m

Skin surface



VS

Ampoule + Hands

0~15 μ m

Skin surface



Ampoule
(Shown in red)


VS

Ampoule + Serum Stick (EP)

35~50 μ m

Up to **3x deeper** absorption
when using the device

Skin surface



Proven deep
skin penetration
of ampoule
(shown in red)

Skin Penetration Test by Human Skin Clinical Trial Center Co., Ltd.

SERUM STICK

A hygienic and convenient stick-type device
for powerful skincare

Portability

Designed in a compact stick form, it's easy to carry anywhere.
Enjoy quick and effortless skincare even on the go.

Hygienic Use

The stick format avoids direct hand contact,
keeping it hygienic and mess-free, no drips, no sticky residue.

Convenience

With just one simple use a day, you can experience effective skin
benefits — perfect for those with busy lifestyles.

High-Concentration Serum Delivery

Powered by EP (Electroporation) technology,
it delivers high-concentration serum deep into your skin—
targeting where it's needed most for visible, lasting results.

“ Serum Stick enhances your daily skincare routine
efficient, effective,

with simple use and powerful results
and enjoyable. ”



Perfect care with one drop,
3 types of highly concentrated serums



VIVA SERUM

#SkinTexture #VitaminC #SpotRemoval

Regulates melanin for a clear,
radiant complexion.

Vitamin C peptides and Bakuchiol
improve skin tone, reduce dark spots,
and brighten the skin.



CENTELLIEVE SERUM

#SkinBarrierStrengthening #SkinSoothing

Customized care for sensitive skin
that strengthens the skin barrier.

Centella Asiatica and Madecassoside
soothe and repair damaged skin,
restoring its health.



STEMPID SERUM

#Elasticity #StemCells #Regeneration

Promotes collagen production
and revitalizes skin cells.

Stem cells regenerate skin,
while peptides firm
for youthful elasticity.